

MPR Upcurve Basket Forceps

10025



Description

Forceps for excising tissue when resecting a damaged meniscus or sculpting remaining meniscal rim into the desired shape. Features a special “basket” tip to ensure a secure grip on tissue and a long, narrow shaft to maximize precision and reduce obstruction, with a slight curve to improve access under the femoral condyles and to the posterior compartment of the knee.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and designed for optimal performance with minimal risk of trauma.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	9" (22.86cm)
Tip Width	0.16" (4mm)
Tip Spread	0.49" (12.5mm)
Hole Length	0.28" (7mm)
Hole Width	0.12" (3mm)
Bend	15°

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.