



Carpal Tunnel Retractor

Description

A retractor used during carpal tunnel-related procedures after the incision has been made, to elevate tissue off of the transverse ligament, allowing surgery to be performed under direct visualization. Features a slightly hooked, stepped tip to help raise the proximal tissue, and a curved handle end for reduced fatigue and enhanced precision.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	5.06" (25.4cm)
Tip Length	0.63" (1.59cm)
Point Length	0.13" (0.32cm)
Point Angle	90°

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

