



MPR Knee Cannula w/ Holes

10766



Description

A cannula developed for use in knee arthroscopic procedures. Features a knurled knob for secure, stable handling and easy adjustment, and a pair of holes near the tip to allow inflow to escape if the end becomes lodged in tissue or closed-off due to valgus or varus stress. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	3.77" (9.58cm)
Shaft Diameter	0.22" (0.56cm)
Hole Diameter	0.2" (0.51cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.