

MPR Wagner Mallet

Description

An orthopedic mallet, typically used with a chisel to cut and shape bone. Features a ridged handle for a secure, comfortable grip, and interchangeable steel and plastic caps. Made from solid stainless steel construction for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	9.75" (24.8cm)
Overall Weight	16.93oz (480g)
Head Diameter	1.46" (3.7cm)
Head Weight	10.58 (300g)

Accessories

21178	Pair of Metal Caps
21179	Pair of Plastic Caps
21180	Wrench for Changing Caps



For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.