



Cannulated Screw Driver

Description

A driver used in conjunction with Instrument Makar mini screws, with a cannulated shaft. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	5.75" (14.61cm)
Shaft Diameter	0.24" (0.61cm)
Tip Diameter	0.1" (0.25cm)
Hole Diameter	0.43" (1.1cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

