



## MPR Slotted Mallet, Metal Handle

MPR-1120-1000



### Description

A mallet with a slotted head for inserting and removing reamers. Features a fenestrated handle for reduced weight and fatigue from use, and reduced heat retention after sterilization. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

|                |                 |
|----------------|-----------------|
| Overall Length | 8.75" (22.23cm) |
| Overall Weight | 37.6oz (1,066g) |
| Handle Length  | 4.49" (11.4cm)  |
| Head Diameter  | 1.75" (4.45cm)  |
| Slot Width     | 0.56" (1.43cm)  |
| Slot Depth     | 0.84" (2.13cm)  |

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).