



MPR ACL Reconstruction Forked Tap

MPR-12001



Description

A forked tap developed for proper placement and insertion of the Biologically Quiet™ Staple during ACL reconstruction procedures, used to prepare the implant site by impacting it into the previously compacted ceiling of the femoral socket. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	9" (22.86cm)
Grip Diameter	1" (2.54cm)
Shaft Diameter	0.28" (0.71cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.