



MPR ACL Reconstruction Obturator

MPR-12004



Description

An obturator used with the MPR ACL Reconstruction Cannula (12003), for placement and insertion of the Biologically Quiet™ Staple during ACL reconstruction procedures. Features a knurled gripping end for secure, stable handling. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	4.13" (10.48cm)
Shaft Diameter	0.39" (0.99cm)
Tip Length	0.29" (0.74cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.