



MPR K-Wire Tractors



MPR-121-02

MPR-121-03

Description

Also known as traction bows, Kirschner bows, or Kirschner extension bows. Instruments for holding K-Wires in tension during orthopedic procedures. The serrated clamps at the ends of the arms ensure a firm grip on the wire, and the handles feature five holes for attaching hooks for suspension. MPR Orthopedics offers K-Wire Tractors in two different sizes; both are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization, and come with a set of 3 hooks.

Use with threaded K-Wires to ensure a proper grip.

Product Specifications

#	LENGTH	ARM LENGTH	HANDLE WIDTH	MAX SPREAD
MPR-121-02	10.5" (26.67cm)	7.48" (19cm)	2.94" (7.47cm)	9.84" (25cm)
MPR-121-03	12.75" (32.39cm)	9.45" (24cm)	2.94" (7.47cm)	11.75" (29.85cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.