



MPR Moreland Hook

MPR-6020-12



Description

A hook used during revision procedures for cement removal/extraction, particularly removal of cement restrictors from the femur. Features a T-handle to ensure a secure, comfortable grip and optimal performance, a long, thin shaft with a tapered end for deep access and precision, and a curved hooked tip. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA. Designed by Dr. John R. Moreland, M.D.

Product Specifications

Overall Length	17" (43.18cm)
Handle Width	3.15" (7.99cm)
Shaft Length	13.25" (33.66cm)
Shaft Diameter	0.23" (0.58cm)
Tip Length	0.27" (0.69cm)
Tip Width	0.25" (0.64cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.