



MPR Straight Ball Spike

MPR-7117



Description

A straight ball spike for reducing and holding fractured bones in place. The sharp point easily digs into the bone tissue, eliminating the risk of slippage, and the faceted Radel®* handle prevents rolling when placed down. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	11.81" (30cm)
Handle Length	5.51" (14cm)
Handle Diameter	1" (2.54cm)
Shaft Length	5.1" (12.95cm)
Shaft Diameter	0.31" (0.8cm)
Ball Diameter	0.31" (0.79cm)
Tip Length	0.1" (0.25cm)
Tip Diameter†	0.1" (0.23cm)

†measured at base

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com