



MPR Femoral Head Driver

MPR-9027-58



Description

An impactor for driving a femoral implant into the socket during Total Hip Arthroplasty procedures. Features an ergonomically designed aluminum handle for a secure, comfortable grip and easy cleaning/sterilization, a striking cap, and a removable Delrin® head, with a depressed center to ensure proper engagement with the implant.

Made in USA.

Product Specifications

Overall Length	6.63" (16.19cm)
Handle Length	4" (10.16cm)
Handle Diameter*	0.63" (16mm)
Head Length	1" (1.54cm)
Head Diameter	1.5" (3.81cm)

*measured at widest point, also applies to striking cap

Accessories

MPR-9027-58H	MPR Replacement Head for Femoral Driver
--------------	---

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.