



MPR Adjustable Leg Rest

MPR-9390



Description

A standalone leg rest, especially useful in casting, amputation, and/or skin prepping. Its adjustable height makes it highly versatile, and the wide square base ensures stability during use. Made from stainless steel for durability, longevity, and easy cleaning/disinfecting.

Made in USA.

Product Specifications

Min. Height	7" (17.78cm)
Max. Height	10" (25.4cm)

Cleaning & Maintenance

Between uses, wipe device down with a hospital- and/or EPA-approved germicide, or a **10:1** bleach solution or quaternary ammonia solution, then rinse and let dry. *Do not* use phenolic cleaners or abrasive materials.

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.