

MPR Clark Sockets



MPR-CS



MPR-SRC



MPR-CS-AIO

Description

Clark-style sockets for mounting positioning devices such as leg holders and hip positioners to surgical tables. All are designed to fit any North American surgical table rail, and accept standard $\frac{5}{8}$ " round posts and round bar accessories up to $1\frac{1}{4}$ ". Additionally, the MPR All-in-One Clark Socket (MPR-CS-AIO) can accept blade-style posts up to 1" wide, and fits onto Euro-size rails as well.

Two different attachment styles are available; the MPR All-in-One Clark Socket and MPR Stainless Steel Clark Socket (MPR-CS) attach via sliding onto the rail, whereas the MPR Over Side Rail Clark Socket (MPR-SRC) can be slipped directly onto any point on the rail.

All are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization. Made in USA.

Product Specifications

#	DESCRIPTION	LENGTH	HANDLE LENGTH	WIDTH
MPR-CS	MPR Stainless Steel Clark Socket	3.88" (9.84cm)	5.75" (14.61cm)	2.81" (7.14cm)
MPR-SRC	MPR Over Side Rail Clark Socket	4" (10.16cm)	4.5" (11.43cm)	3.5" (8.89cm)
MPR-CS-AIO	MPR All-in-One Clark Socket	4.75" (12.07cm)	5.63" (14.29cm)	3.63" (9.21cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.