



MPR Instrument Tables, Height-Adjustable



Description

Versatile instrument tables for use in operating rooms. These tables are designed to straddle the OR table, and feature two locking casters and two swivel casters, all with a 3" heavy duty design, and a removable crank handle which can be attached to adjust the height from either side. MPR Orthopedics offers five different models in a variety of sizes. All are made from fully welded 16-gauge stainless steel and have a height range of 36"–56" (91.44cm–142.24cm).

Made in USA.

Product Specifications

#	LENGTH	WIDTH	LEG CLEARANCE
MPR203628	36" (91.44cm)	20" (50.8cm)	28" (71.12cm)
MPR243628	36" (91.44cm)	24" (60.96cm)	28" (71.12cm)
MPR244840	48" (121.92cm)	24" (60.96cm)	40" (101.4cm)
MPR246052	60" (152.4cm)	24" (60.96cm)	52" (132.08cm)
MPR463426	34" (86.36cm)	46" (116.84cm)	26" (66.04cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.