



MPR Surgical Hand Table

MPR-P740



Description

A small table for providing additional support for hand surgeries and other procedures requiring a stable platform. The pointed legs ensure steady placement by digging into the table pad, and its lightweight, stackable construction simplifies setup and storage.

Made from solid stainless steel for longevity, durability, and easy cleaning/disinfecting.

Made in USA.

Product Specifications

Length	8" (20.32cm)
Width	6.25" (15.88cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.