



MPR Vise Grips/Pliers

Description

Pliers used in a number of procedures, such as revision, extraction, and explanting. Their flat designs are optimized for comfortable one-handed operation with a minimal profile, and feature a variety of jaw styles to suit your needs/preferences, and up to two threaded inserts for attaching a slap hammer, such as the MPR Slap Hammer (MPR831). All are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

MPR Curved Jaw Vise Grips/Pliers



#	LENGTH	JAW LENGTH	BACK INSERT	SIDE INSERT
MPR820	5" (12.7cm)	0.88" (2.22cm)	Yes	Yes
MPR822	8" (20.32cm)	1.38" (3.49cm)	Yes	Yes
MPR823	10" (25.4cm)	1.5" (3.81cm)	Yes	Yes

MPR Straight Jaw Vise Grips/Pliers



#	LENGTH	JAW LENGTH	BACK INSERT	SIDE INSERT
MPR819	10" (25.4cm)	1.38" (3.49cm)	No	No
MPR824	7" (17.78cm)	1.38" (3.49cm)	Yes	Yes
MPR825	10" (25.4cm)	1.38" (3.49cm)	Yes	Yes

MPR Long Nose Vise Grips/Pliers



#	LENGTH	JAW LENGTH	BACK INSERT	SIDE INSERT
MPR826	6" (15.24cm)	2.5" (5.72cm)	Yes	No
MPR827	9" (22.86)	2.75" (6.99cm)	Yes	Yes
MPR828	11" (27.94cm)	2.5" (5.72cm)	Yes	Yes

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.