



MPR Double-Handed Bone Curettes



MPR00403301700

MPR00403301800

MPR00403301900

Description

Curettes for debriding and excising bone tissue during revision procedures. Their extended designs allow them to be held with two hands for maximal stability and effectiveness, and feature knurled grips for secure, comfortable handling. MPR Orthopedics offers three models with differing tip sizes. All are made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

#	SIZE	LENGTH	GRIP LENGTH	TIP LENGTH	TIP WIDTH*
MPR00403301700	Small	15" (38.1cm)	4.31" (10.95cm)	0.75" (1.91cm)	0.31" (0.79cm)
MPR00403301800	Medium	15" (38.1cm)	4.56" (11.59cm)	0.88" (2.22cm)	0.5" (1.27cm)
MPR00403301900	Large	15.13" (38.42cm)	4.25" (10.8cm)	1.13" (2.86cm)	0.63" (1.59cm)

*measured at widest point

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.