



MPR Width Plates



MPR00880200100



MPR00880200200



MPR00880200300



MPR00880200400

Description

Width plates designed for use with Zimmer* dermatomes that use a two-screw attachment style, to determine the width of skin grafts. MPR Orthopedics offers four different sizes for a variety of width options. Made from high-grade aluminum for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

#	LENGTH	OVERALL WIDTH	THICKNESS	GRAFT WIDTH
MPR00880200100	2.13" (5.4cm)	4.75" (12.07cm)	0.13" (0.32cm)	1" (2.54cm)
MPR00880200200	2.13" (5.4cm)	4.75" (12.07cm)	0.13" (0.32cm)	2" (5.08cm)
MPR00880200300	2.13" (5.4cm)	4.75" (12.07cm)	0.13" (0.32cm)	3" (7.62cm)
MPR00880200400	2.13" (5.4cm)	4.75" (12.07cm)	0.13" (0.32cm)	4" (10.16cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.