



Stand-Alone Leg Prepper

Description

A standalone leg rest, usable with any table, gurney, or stretcher. Features an arm and head that can be easily raised, lowered, and rotated as needed, and a base designed for optimal stability. Made from stainless steel for durability, longevity, and easy cleaning/disinfecting.

Made in USA. Some assembly may be required.

Product Specifications

Min Height	44" (111.76cm)
Max Height	72" (182.88cm)
Width	19" (48.26cm)
Depth	19" (48.26cm)
Weight Limit	125lbs (56.7kg) (approximated)



Accessories

#	DESCRIPTION
MPR054-P	MPR Reusable Replacement Pad

Cleaning & Maintenance

Between uses, wipe device down with a hospital- and/or EPA-approved germicide, or a **10:1** bleach solution or quaternary ammonia solution, then rinse and let dry. Avoid using phenolic cleaners or abrasive materials.

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.