



MPR Browne Deltoid Retractors



MPR1670-01

MPR1670-02

Description

Retractors used during shoulder orthopedic and revision procedures to retract the deltoid muscle and surrounding tissue. Their ankh-shaped handles feature a ring which can be used to attach weights or as a finger-hole, and in conjunction with their straight shafts, ensure optimal handling and leverage, while the concave, sharply angled blades easily slide under muscle tissue and engage directly with the bone. MPR Orthopedics offers two different sizes, to suit varied patient requirements and your preferences.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

#	DESC.	LENGTH	HANDLE WIDTH	BLADE LENGTH	BLADE WIDTH
MPR1670-01	Small	10.5" (26.67cm)	2.75" (6.99cm)	2.75" (6.99cm)	1.63" (4.13cm)
MPR1670-02	Large	10.75" (27.31cm)	2.75" (6.99cm)	2.75" (6.99cm)	2.13" (5.4cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.