



MPR Moreland Poly Implant Extractor Sleeve

MPR2420-05S



Description

A replacement sleeve for DePuy's* Poly Implant Extractor Handle (2420-05-000). Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. John R. Moreland, M.D.

Product Specifications

Overall Length	5.5" (13.97cm)
Grip Length	3.5" (8.89cm)
Grip Diameter	0.75" (1.91cm)
End Diameter	1.5" (3.81cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.



*DePuy is a registered trademark of that corporation.