

Hip Retractor, 5" Handle

MPR3319-05



Description

A retractor specialized for use in hip surgery. The ergonomic silicone handle is designed to ensure a secure, comfortable grip with minimal risk of slipping, and the curved blade provides optimal access, with a thick hooked tip for solid tissue engagement.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	13" (33.02cm)
Handle Length	5" (12.7cm)
Blade Length	6" (15.24cm)
Blade Width	1" (2.54cm)
Tip Length	1" (2.54cm)
Tip Width*	0.25" (0.64cm)

*measured at base

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.