



MPR 3mm Probe Hook

MPR3856



Description

A hook with a blunt, teardrop-shaped tip for tissue probing and minute tissue manipulation and/or retraction during general surgical procedures. The handle is knurled to ensure a secure, comfortable grip with minimal risk of fatigue from prolonged use, and the slender shaft provides maximal access as well as delicate control and tactile feedback. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	7.25" (18.42cm)
Handle Length	3.75" (9.53cm)
Handle Diameter	0.19" (0.47cm)
Hook Length	0.12" (0.12cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.