



MPR Flag Splitter

Description

A splitter primarily used during revision procedures to fragment the remaining cement for removal after a femoral component has been removed. Features a tip with a sharp angled blade, a blunt guiding point to help navigate the tip along the cortices to prevent perforation of the femur, and an ergonomically designed Radel®* handle with a strike plate and tip alignment mark.

Made in USA.

Product Specifications

Overall Length	17.75" (45.09cm)
Handle Length†	5.62" (14.28cm)
Handle Diameter	1.3" (3.29cm)
Shaft Length	11.13" (28.28cm)
Shaft Diameter	0.25" (0.64cm)
Tip Length	1.01" (2.57cm)
Tip Width	0.31" (0.8cm)
Blade Length	0.87" (2.22cm)
Blade Width	0.12" (0.3cm)
Blade Angle	58°

†including strike plate

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com

MPR906

