



MPR Acetabular Component Gripper

Description

A self-retaining gripper used in conjunction with a slap hammer to remove a cup with an intact polyethylene liner. Comprises a straight and a curved arm, both with outward-facing points, and a speed-lock mechanism. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	12.66" (32.15cm)
Grip Length	3.81" (9.68cm)
Grip Diameter	0.62" (1.58cm)
Jaw Length	2.21" (5.62cm)
Point Length	0.28" (0.72cm)

Product Compatibility

#	DESCRIPTION
MPR831	MPR Slap Hammer
MPR939	MPR Small Slap Hammer

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

