



MPR Acetabular Component Forceps

Description

A specialized forcep used to extract a polyethylene acetabular component liner during a liner exchange. Comprises a straight arm with a single-pronged hook tip and a curved arm with a forked tip, ensuring optimal access and gripping ability. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	15" (38.1cm)
Handle Length	11.63" (29.55cm)
Handle Width	0.39" (0.52cm)
Straight Jaw Length	3.37" (8.55cm)
Curved Jaw Length	3.48" (8.85cm)
Max. Opening Angle	52°

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

