

## MPR Closed Loop Extractor

MPR946



### Description

A loop extractor for removing a one-piece prosthesis stem during revision and extraction procedures. The loop design fits neatly around the the head of the prosthesis and the angled tip provides leverage for the extraction process. Attaches to a slap hammer via the quick connector on the shaft. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

|                |                 |
|----------------|-----------------|
| Overall Length | 5.87" (14.92cm) |
| Overall Width  | 2.17" (5.52cm)  |
| Tip Length     | 1.9" (4.84cm)   |
| Tip Width      | 1.04" (2.63cm)  |
| Tip Angle      | 30°             |
| Wire Diameter  | 0.23" (0.59cm)  |

### Product Compatibility

| #      | DESCRIPTION                       |
|--------|-----------------------------------|
| MPR950 | MPR Femoral Extractor Slap Hammer |

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).