



MPR Lucae Bone Mallet, Stainless Steel

21120



Description

A bone mallet for use during osteotomy and other revision procedures. Features an ergonomically designed handle for a secure comfortable grip, convex striking faces, and solid stainless steel construction for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	7.5" (19.05cm)
Overall Weight	7oz (200g)
Head Diameter	0.75" (1.9cm)
Head Weight	3oz (80g)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.