



MPR Narrow-Tip Mallet

Description

A specialized orthopedic mallet ideal for hand procedures or other areas with small structures. Features an ergonomic phenolic resin handle for a secure, comfortable grip, and a head with one square striking face and one protruding narrow circular striking face, for use in delicate areas, or those with limited access. Made from stainless steel for durability, longevity, and easy cleaning/sterilization. The handle can be safely steam-sterilized up to 250°F/121°C.



Product Specifications

Overall Length	8" (20.32cm)
Overall Weight	7oz (198.45g)
Square Face Width	0.71" (1.8cm)
Tip Diameter	0.28" (0.7cm)
Head Weight	3oz (85.05g)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.