



MPR ACL Reconstruction Socket Compactor

MPR-12005



Description

An instrument used to compact the ceiling of the femoral socket prior to insertion of the Biologically Quiet™ Staple during ACL reconstruction procedures, to ensure solid fixation. Features a T-handle for secure, stable handling. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

Product Specifications

Overall Length	8.48" (21.54cm)
Handle Width	2" (5.08cm)
Shaft Diameter	0.25" (0.54cm)
Tip Length	0.56" (1.43cm)
Tip Diameter	0.43" (1.1cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.