



MPR Femoral Head Extractor

MPR-1440-1010



Description

A corkscrew-type extractor for removing the femoral head during revision procedures. The sharp, threaded tip is designed to provide a secure grip and cleanly drill into the femoral head. Easily attaches to the MPR T-Handle (MPR980), or any other handle with a Zimmer/Hall-style connector; verify fit before use.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	9.15" (23.24cm)
Shaft Length	6.83" (17.34cm)
Shaft Diameter	0.31" (0.79cm)
Tip Length	1.33" (3.37cm)
Tip Diameter*	0.4" (1.02cm)

*including threads

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.