



MPR Extra Leverage Femoral Neck Elevator

MPR-2598-765



Description

An instrument used during hip procedures such as Total Hip Arthroplasty to elevate and expose the femoral head. The handle is knurled to ensure a secure, comfortable grip, and can accommodate either one or both hands, and the wide, flat blade is bowed slightly to provide a greater amount of leverage than other femoral neck elevators, with a serrated tip guaranteeing firm engagement with tissue.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	17.87" (45.4cm)
Handle Length	8.86" (22.5cm)
Handle Diameter	0.99" (2.46cm)
Blade Length	9.02" (22.9cm)
Blade Width	1.5" (3.8cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.