

MPR Femoral Neck Elevator

MPR-7803-74



Description

An instrument used during hip procedures such as Total Hip Arthroplasty to elevate the proximal femur and expose the femoral head. The bent handle and blade grant the user a high degree of leverage, while the thin profile enhances access and the two-pronged tip ensures a secure grip on the femur. Designed for optimal performance with minimal risk of trauma.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	15.97" (40.56cm)
Handle Length	3.42" (8.69cm)
Handle Width	1.3" (3.3cm)
Blade Length	4.7" (11.94cm)
Blade Width	1.58" (4cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.