



## MPR 13" Narrow Tip Retractor

MPR1053



### Description

A retractor primarily used during hip procedures to retract soft tissue and provide greater exposure of the joint. Its long, narrow blade/tip is designed for deep access to retract a large amount of tissue with minimal visual and physical obstruction. Features a hole in the handle for attaching weights.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

### Product Specifications

|                |               |
|----------------|---------------|
| Overall Length | 13" (33.02cm) |
| Tip Width      | 0.47" (1.2cm) |
| Point Width    | 0.2" (0.52cm) |

### For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at [www.mprortho.com](http://www.mprortho.com).