



MPR Replacement Clamp for MPR330



Description

A replacement clamp for the MPR330 MPR Carpal Arm and Hand Table. Designed to fit any American OR table rail, its lightweight aluminum construction ensures durability while being easy to install onto the hand table and providing sturdy, consistent alignment and support.

Made in USA.

Product Specifications

Clamp Width	1.57" (4cm)
Clamp Height	0.87" (2.2cm)
Jaw Opening	0.47" (1.2cm)
Mounting Slot Width	0.24" (0.6cm)

Cleaning & Sterilization

Wipe down thoroughly with a hospital- and/or EPA-approved cleaner and disinfectant prepared according to the manufacturer's instructions. If using a cleaning spray, apply the spray to a cloth and wipe instead of spraying the device directly. Avoid phenolic or bleach-based cleaning products, or any harsh abrasives. Use another clean, dry cloth to remove the cleaning fluid, making sure that there is no pooling, then let dry.

DO NOT use an autoclave to sterilize this device.

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.