



Stainless Steel Anesthesia Cart

MPR523



Description

An anesthesia cart with a large tray, and four drawers and lower shelf for extra storage. Made from heavy-gauge solid stainless steel for durability and reliability, and two locking casters.

Made in USA.

Product Specifications

Length	16" (40.64cm)
Width	20" (50.8cm)
Height	34" (86.36cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.