



MPR Moreland Slotted Mallets



MPR928

MPR938

Description

Orthopedic mallets with slotted heads for inserting and removing reamers. The handles are made from Radel^{®†} and are ergonomically designed for a comfortable, secure grip, and feature stainless steel strike plates at the ends.

MPR Orthopedics offers Moreland Slotted Mallets in two different sizes. The larger mallet features a removable Delrin[®] cap on one of the striking faces, offering greater versatility in use.

Made in USA. Original design by Dr. John R. Moreland, M.D.

Product Specifications

#	DESCRIPTION	LENGTH	HEAD WIDTH	HEAD DIAMETER	WEIGHT
MPR928	Moreland Slotted Mallet	9" (22.86cm)	2.75" (7cm)*	1.56" (40mm)	1.6lb (726g)
MPR938	Moreland Small Slotted Mallet	8.13" (20.64cm)	2.5" (6.35cm)	1.31" (33.25mm)	1.4lb (635g)

*with cap

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

[†]Compared with other materials commonly used on handles for surgical instruments, Radel[®] is engineered for harsh environments, and can be cleaned, disinfected, and steam sterilized over 1,000 cycles without significant loss of properties.