



MPR Femoral Canal Finder



Description

An instrument used in Total Hip Arthroplasty procedures using a direct anterior approach to locate and prepare the femoral canal for stem broaching. Its tapered, slightly bent design is optimized for easy entry into the canal, with a smooth working end to minimize undesired tissue damage. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	9.69" (24.6cm)
Handle Width	4.92" (12.5cm)
Handle Diameter	0.59" (1.5cm)
Bend	27°
Tip Diameter	0.16" (0.4cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.