# MPR ACL Reconstruction Cannula

12003



## Description

A cannula developed for proper placement and insertion of the Biologically Quiet™ Staple during ACL reconstruction procedures. Features a knurled gripping end for secure, stable handling. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Used in conjunction with the corresponding obturator (MPR-12004).

Made in USA.

Designed by Dr. Lanny Johnson, M.D.

## Product Specifications

| | |
|---|---|
| Overall Length | 3.44” (8.73cm) |
| Shaft Diameter | 0.44” (1.11cm) |
| Hole Diameter | 0.4” (1cm) |

## For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com