MPR Flexible Osteotomes, Bayonet Style



Description

Flexible blades used in conjunction with a compatible handle for osteotomies and implant loosening. All are made from solid stainless steel for maximal effectiveness, and feature a bayonet-style connector.

Made in USA. Available either sterile or nonsterile. Intended for single use *only*; safely dispose after use.

Product Specifications

#	BLADE LENGTH	BLADE WIDTH
MPR954	3” (7.62cm)	0.31” (0.8cm)
MPR955	3” (7.62cm)	0.39” (1cm)
MPR956	3” (7.62cm)	0.47” (1.2cm)
MPR957	3” (7.62cm)	0.79” (2cm)

#	BLADE LENGTH	BLADE WIDTH
MPR960	5” (12.7cm)	0.31” (0.8cm)
MPR961	5” (12.7cm)	0.39” (1cm)
MPR958*	3” (7.62cm)	0.47” (1.2cm)
MPR959*	5” (12.7cm)	0.79” (2cm)

*curved end

Product Compatibility

#	DESCRIPTION
MPR948	MPR Flexible Osteotome Handle, Bayonet Style

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com