
MPR Reamer Hole Gauge

MPR210599000



Description

A measurement tool for checking and verifying reamer size, detecting wear or undersizing on reamers, and ensuring accurate canal sizing when fitting impants. The durable stainless steel construction minimizes wear on the gauge, guaranteeing consistently reliable measurements with years of repeated use, and is easy to clean.

Made in USA.

Product Specifications

Length	11” (27.94cm)
Width	4.75” (12.07cm)
Minimum Hole Ø	0.31” (8mm)
Maximum Hole Ø	0.89” (22.5mm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.