MPR Sharp Hook

MPR-319.39





Description

A versatile hook with a sharp tip, used during a number of procedures including revision and fracture repair operations, for manipulating delicate bone fragments, checking alignments, and removing tissue from the heads and threads of screws. The hook is designed for maximal dexterity and precision with minimal obstruction, and the octagonal handle is ridged for a secure, comfortable grip. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	6” (15.24cm)
Handle Length	4.6” (11.68cm)
Handle Diameter	0.22” (0.55cm)
Hook Length	0.83” (2.12cm)
Hook Diameter*	0.05” (0.13cm)
Tip Width	0.18” (0.46cm)

*averaged

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com