MPR Cartilage Scissors

MPR247700000



Description

Specialized scissors for cutting cartilaginous tissue. Their curved blades provide easier access to anatomically challenging areas and allow for precise cutting of small amounts of tissue, and have serrated edges to maximize their cutting ability and eliminate risk of tissue slippage. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	8” (20.32cm)
Blade Length	2” (5.08cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com