MPR Moreland Straight Gouge

MPR910



Description

A gouge primarily used for cement removal and extraction during revision procedures, with an ergonomic Radel®* handle for a secure, comfortable grip, and a strike plate for use with a mallet, featuring an engraved alignment mark to guide tip orientation. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA. Original design by Dr. John R. Moreland, M.D.

Product Specifications

Overall Length	17.25” (43.82cm)
Handle Length†	5.62” (14.28cm)
Handle Diameter	1.59” (4.04cm)
Shaft Length	12.25” (31.12cm)
Shaft Diameter	0.25” (0.63cm)
Tip Length	1” (2.54cm)
Tip Width	0.38” (0.95cm)

†including strike plate

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

*Compared with other materials commonly used on handles for surgical instruments, Radel® is engineered for harsh environments, and can be cleaned, disinfected, and steam-sterilized over 1,000 cycles without significant loss of properties.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com