MPR Lead-Filled Mallet

21124



Description

A lead-filled and -faced mallet, designed to reduce damage from strikes. Its ergonomically designed handle ensures a secure, comfortable grip, and the stainless-steel-clad head protects the soft lead interior, ensuring its durability and longevity.

This mallet should only be used to strike flat surfaces, as hitting sharp objects can cause lead to flake off.

Made in USA.

Product Specifications

Overall Length	7.75″ (19.5cm)
Overall Weight	7oz (200g)
Head Diameter	0.87″ (2.2cm)
Head Weight	5oz (150g)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com