MPR Microfracture Mallet

7777-02



Description

A lightweight orthopedic mallet, ideal for microfracture repair procedures. The design features a ridged handle for a secure grip, and is made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	9” (22.86cm)
Handle Length	4.63” (12.6cm)
Head Diameter	1.13” (2.8cm)
Weight	13oz (363g)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com