MPR Stainless Steel Spoon

MPR-5900-N



Description

A spoon used as a curette for removal of cement fragments and other material during revision procedures. Its handle is designed to be comfortable to hold and reduce fatigue from use, and the broad, egg-shaped scoop allows for removal of both large and minute amounts of material, and offers a high degree of precision. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

Overall Length	7.28” (18.5cm)
Handle Length	4.72” (12cm)
Scoop Length	1.65” (4.2cm)
Scoop Width	2.36” (6cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com