MPR “Thor’s Hammer” Mallet

MPR-1163-1000



Description

An orthopedic mallet designed to resemble Thor’s hammer Mjölnir from the Marvel Studios films. Features a steam-sterilizable propylux handle with a stainless steel core, ensuring a secure, comfortable grip with minimal risk of slippage. Made from stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	8.88″ (22.54cm)
Overall Weight	35.27oz (1,000g)
Handle Length	7” (17.78cm)
Head Length	1.44” (4.23cm)
Head Width	3.2 (8.13cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com