MPR Blunt Narrow Cobra Retractor

MPR-1440-1140



Description

A retractor used during hip procedures utilizing an anterior approach to assist in exposing the joint. Features a curved blade to provide deep tissue access, and a bent tip to help ensure a secure grip on the tissue. Features four holes in the handle for attaching weights.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

Overall Length	12” (30.48cm)
Blade Length	4” (10.16cm)
Blade Width	0.89” (2.27cm)
Tip Length	1.2” (3.05cm)
Tip Width	0.48” (1.21cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com